

Kodak DIRECTVIEW Total Quality Tool

Overview

Xiaohui Wang, PhD
Richard Van Metter, PhD
Eastman Kodak Company
Rochester, NY

HEALTH IMAGING
A BETTER VIEW OF LIFE.



KODAK Digital Capture Solutions

Milestones



Manual Acceptance Testing and Quality Control

Technical and Scientific Bulletin

Guidelines for Acceptance Testing and Quality Control

Kodak DirectView CR 800 System
and
Kodak DirectView CR 900 System

- Preliminary testing (x-ray machine, displays, printers)
- Inventory & Inspection
- Throughput
- Linearity
- Uniformity & Artifacts
- Erase Function
- Geometry
- Cassette testing (Exposure Response, Uniformity and Artifacts)

Acceptance Testing:

Quality Control Testing:

CT800 Acceptance T1

Yearly	<ul style="list-style-type: none">• Complete acceptance test
Twice yearly	<ul style="list-style-type: none">• Cassette Exposure Response, Uniformity and Artifacts - Test the 10 mR system response and visually check the resulting image from each cassette for uniformity and artifacts.
Monthly	<ul style="list-style-type: none">• Visually inspect all screens for dust and scratches
Weekly	<ul style="list-style-type: none">• Erase all unused cassettes• Verify luminance calibration of workstation displays
Daily	<ul style="list-style-type: none">• Verify printer-processor density calibration



Kodak TQT – CR Testing Matrix

Manufacturing \Leftrightarrow Service \Leftrightarrow User

Exposure Response

- *Linearity & Noise*

Spatial Resolution - MTF

- *Slow scan & fast scan*

- *50% and 95% f_{Nyquist}*

Geometric Accuracy

- *Pixel Spacing, Aspect Ratio, Scan Linearity*

Field Uniformity

Erase

Artifacts:

- *Streaks, Pixel-Position Error, Line-Position Error*


System Electronic Noise




Acceptance Testing and Quality Control

Kodak Total Quality Tool

- **Phantom** - *targets for quantitative analyses*
- **Procedure** - *acquire phantom and flat-field images using controlled exposures*
- **Analysis Software** - *automatic image analyses and decision making*
- **Documentation**
 - **Theory of Operation**
 - **User Guide Addendum**
 - **Quick Users**





Publication No. XXXXXX Rev 01 © Eastman Kodak Company, 2001

Quick Reference Guide for the Kodak DirectView Total Quality Tool

CR SYSTEM TESTING




Starting the Total Quality Tool:
From the CR System main menu, touch
Key Operator, then touch Total Quality Tool.

Closing the Total Quality Tool:


Performing CR System Tests:
DO NOT load the cassette before starting the test.

1. At Total Quality tool main page, touch the button for the test you wish to perform.
2. When you see the "Load cassette..." message, load the appropriate test cassette for the selected test.


Kodak TQT Phantom



18 x 24 cm




24 x 30 cm




35 x 43 cm

Kodak TQT Procedure




0.5 mm Cu
1.0 mm Al

10.0 ± 0.2 mR
@ 80 kVp



≥ 180 cm



Publication No. XXXXXX Rev 01 © Eastman Kodak Company, 2001

Quick Reference Guide for the Kodak DirectView Total Quality Tool

CR SYSTEM TESTING

Starting the Total Quality Tool:
From the CR System main menu, touch **Key Operator**, then touch **Total Quality Tool**.

Closing the Total Quality Tool:
Touch **Main Menu**.

Workflow Tips for CR System Testing:

1. Acquire the Phantom test image and Flat Field test image on separate cassettes.
2. Perform the Phantom Image Test and then the Flat Field Image Test.
3. Use the erased cassette from the Flat Field Image Test to perform the Erase Image Test and the System Noise Test.

Acquiring Test Images:

0.5 mm Cu
1.0 mm Al

10.0 ± 0.1 mR
@ 80 kVp

Lead Apron

Acquiring a Phantom Test Image

1. For the Phantom image, place a cassette in the phantom tray and place the Phantom Test Plate inside the tray.
2. Position the tray so the image can be acquired.
3. Prepare a lead apron.
4. Limit the exposure level for test images to 10.0 ± 0.1 mR @ 80 kVp.
5. Wait 15 minutes between exposure and screen reading.

Performing CR System Tests:
DO NOT load the cassette before starting the test.

1. At Total Quality tool main page, touch the button for the test you wish to perform.
2. When you see the "Load cassette..." message, load the appropriate test cassette for the selected test. FAIL or N/A results will be displayed if an incorrect cassette is used.
3. When the message changes to "Processing Complete", check the results indicators. See "Troubleshooting" on the reverse side for information on failure (red FAIL), approaching specification limit (amber PASS), or N/A results.

Viewing the Results Graph for a Specific Subtest/Cassette:

1. From the Total Quality Tool main page, touch **Results**.
2. Go to the results page for the test you want to see (Results Page 1 for Phantom Image Test, Results Page 2 for all others).
3. Touch one of the cassette size buttons at the bottom of the page.
4. Touch the appropriate cassette size button to change the cassette type from GP to HR (or vice versa).
5. Touch (pic) to the right of the subtest name.

Accessing the Test Data Summary:


1. Touch **Results** on the Total Quality Tool's main page.
2. Touch **Test Data** to access a summary of the data for (up to) the last 13 tests performed.

Exporting the Test Summary Data:

1. Open the CR System door. Insert a 3 1/2" high-density, blank formatted disk into the disk drive and close the door.
2. Touch **Export** to start processing.



KODAK TQT User Interface



Test Result Details

Phantom Image Test Results

Pixel Size 2.0% max: **1.4%**
Aspect Ratio 3.0% max: **2.7%**
Scan Non-Linearity 2.0% max: **0.95%**
Exposure Latitude 100.0 max: **76.4**
Noise Level (@0.19mR) 13.3 max: **11.9**
Noise Level (@2.11mR) 11.3 max: **10.7**
Noise Level (@10.0mR) 10.3 max: **9.12**
Pixel Placement Error 0.25 max: **0.114**
MTF (short direction @ 50% Nyquist) 34.5% min: **59.5%**
MTF (short direction @ 95% Nyquist) 10.3% min: **21.7%**
MTF (long direction @ 50% Nyquist) 34.5% min: **29.5%**
MTF (long direction @ 50% Nyquist) 10.3% min: **10.9%**

Flat Field Image Test Results

Field Non-Uniformity 100max: **74.0**
Galvo Function 5.0E+06 max: **3.7E+06**
Transport Banding 810.0 max: **715.7**
Random Streak 5 max: **2**

Erased Image Test Results


Residual Signal 100.0 max: **Not Available**

System Noise Test Results


BaseLine System Noise 6.0E-06 max: **Not Available**

Select test result to display graph

MTF (long direction @ 50%)



Previous 8 results of MTF












Test Result Details

DirectView

Results Page 1

directview


Phantom Image Test

Pixel Size Error Fast (%):	0.1	
Pixel Size Error Slow (%):	-0.2	
Aspect Ratio Error Left (%):	-1.5	
Aspect Ratio Error Middle(%):	-0.1	
Aspect Ratio Right Left (%):	-0.8	
Aspect Ratio Error Average (%):	-0.2	
Fast-Scan Speed Error (%):	0.14	
Slow-Scan Speed Error (%):	0.04	
Low-Exposure Response Error (CV):	71.5	
Mid-Exposure Response Error (CV):	-5.6	
High-Exposure Response Error (CV):	-40.4	
Low-Exposure Noise (CV):	16.6	
High-Exposure Noise (CV):	6.6	
High-Exposure Noise (CV):	4.3	
Fast-Scan MTF @ 50% Nyquist (%):	44.1	
Fast-Scan MTF @ 95% Nyquist (%):	17.7	
Slow-Scan MTF @ 50% Nyquist (%):	29.8	
Slow-Scan MTF @ 95% Nyquist (%):	6.0	
Pixel Position RMS (pixels)	0.08	

35 x 43 PHANTOM TEST

Date Tested: 12/18/00

Cassette ID: 999999



MTF DATA

PAGE TWO



TEST DATA



35 x 43



24 x 30 GP



18 x 24 GP



MAIN MENU



BACK





Cassette Quality Assurance

DirectView

Cassette Screen Test

Flat Field Image Test

Field Uniformity: READY

Streaks: READY

Speed Deviation: READY

Previous Session

Date Tested:

Time Tested:

Total # Cassettes Tested:

Total # Tests Performed:

Current Cassette

Cassette ID:

Cassette Size:

Cassette Type:

Total # Actuations:

Append to Previous Session:

Flat Field Image Test

Load Cassette...

Phantom Image Test

Processing Complete

Phantom Image Test

Field Size: PASS

Aspect Ratio: PASS

Scan Linearity: PASS

Exposure Response: PASS

Noise: PASS

MTF: PASS

Flat Field Test: PASS

Flat Field Image Test

Field Uniformity: PASS

Line Position: PASS

Distortion: PASS

Chatter: PASS

Streaks: PASS

Erased Image Test

Erased: READY

System Noise Test

System Noise: READY

RESULTS


MAIN MENU

BACK

TEST DATA

MAIN MENU

BACK



Kodak TQT Summary

Kodak DirectView CR Image Quality Tool is used in production, by service and by users.

- Precise and accurate quality control testing
- Highly reproducible quantitative results
- Detects sub-visible changes in CR image quality performance to initiate timely preventive maintenance
- Avoids hours of tedious and labor-intensive effort with a highly automated procedure
- Full data reporting in Excel format

