



MRI for RT Treatment Planning

Yue Cao, Ph.D.


Departments of Radiation Oncology and Radiology
University of Michigan

No relevant financial relationships

Research is supported in part by P01CA59827,
P01CA85878, RO1CA132834, RO1NS064973

MRI Advantages


- Soft tissue differentiation
- Multiple contrasts



Cao RSNA 2010

Post-Gd T1WI

Nasopharyngeal Cancer



Cao AAPM 2011

Body Sites and Tumors

- Brain tumors
 - Primary and metastatic tumors
- Prostate cancers
 - Delineation of whole prostate gland
 - Localization and Delineation of dominant intra-prostatic lesion
- Cervical cancers
 - Brachy therapy
- Liver tumors
- HN tumors
 - Nasopharyngeal cancer

Cao AAPM 2011

M MRI Advantages

- Soft tissue differentiation
- Multiple contrasts
 - conventional contrasts
 - T1 contrast, T2 (or FLAIR) contrast, Post-Gd T1 contrast
 - Advanced contrasts
 - Susceptibility w (T2*), water and fat separation, cortical bone
 - Molecular, metabolic and functional imaging
 - 1H, 31P and 13C spectroscopy imaging
 - DCE and DSC imaging
 - DW and DT imaging
 - Other contrast agents, e.g., SPIO, Eovist,
 - Hyperpolarized 3He and 13C
- Localization, delineation, and characterization of tumors and normal organs
- Integration of target definition and Tx assessment

Cao AAPM 2010

M Technology Advancements

- High field magnet
- Parallel imaging
- Large Bore size (70 cm)
- Multi-RF transmission
- RF-shimming
- RF coil array/TimCT
- Robust motion suppression pulse sequence



Cao AAPM 2011

M RT Treatment Planning

- Signals, fast acquisition, high resolution, 3D
- RT compatible, embolization equipment
- More uniform RF distribution, e.g., in the liver
- Uniform signal intensity
- Extended coverage and continuous scan like CT
- Better images for motion organ, e.g., liver, HN during swollen



Cao AAPM 2011

M 3D Volumetric T2W image

1x1x1 mm³ resolution on 3T

Cao AAPM 2011

M Technology Advancements

University of Michigan
Medical School

- High field magnet
- Parallel imaging
- Large Bore size (70 cm)
- Multi-RF transmission
- RF-shimming
- RF coil array/TimCT
- Robust motion suppression pulse sequence



Cao AAPM 2011

M RT Treatment Planning

University of Michigan
Medical School

- Signals, fast acquisition, high resolution, 3D
- RT compatible, embolization equipment
- More uniform RF distribution, e.g., in the liver
- Uniform signal intensity
- Extended coverage and continuous scan like CT
- Better images for motion organ, e.g., liver, HN during swollen



Cao AAPM 2011

M Technology Advancements

University of Michigan
Medical School

- High field magnet
- Parallel imaging
- Large Bore size (70 cm)
- Multi-RF transmission
- RF-shimming
- RF coil array/TimCT
- Robust motion suppression pulse sequence



Cao AAPM 2011

M RT Treatment Planning

University of Michigan
Medical School

- Signals, fast acquisition, high resolution, 3D
- RT compatible, embolization equipment
- More uniform RF distribution, e.g., in the liver
- Uniform signal intensity
- Extended coverage and continuous scan like CT
- Better images for motion organ, e.g., liver, HN during swollen



Cao AAPM 2011

M Technology Advancements

University of Michigan
Medical School

- High field magnet
- Parallel imaging
- Large Bore size (70 cm)
- Multi-RF transmission
- RF-shimming
- Modular RF coil arrays/TimCT
- Motion suppression pulse sequence

Cao AAPM 2011



M RT Treatment Planning

University of Michigan
Medical School

- Signals, fast acquisition, high resolution, 3D
- RT compatible, embolization equipment
- More uniform RF distribution, e.g., in the liver
- Uniform signal intensity
- Extended coverage and continuous scan like CT
- Better images for motion organ, e.g., liver, HN during swollen

Cao AAPM 2011



M Technology Advancements

University of Michigan
Medical School

- High field magnet
- Parallel imaging
- Large Bore size (70 cm)
- Multi-RF transmission
- RF-shimming
- RF coil arrays/TimCT
- Robust motion suppression pulse sequences

Cao AAPM 2011



M RT Treatment Planning

University of Michigan
Medical School


- Signals, fast acquisition, high resolution, 3D
- RT compatible, embolization equipment
- More uniform RF distribution, e.g., in the liver
- Uniform signal intensity
- Extended coverage and continuous scan like CT
- Reduced organ movement imaging artifacts

Cao AAPM 2011



M Motion Sensitive?


University of Michigan
Medical School




Cao AAPM 2011

M Motion Suppression

University of Michigan
Medical School



T2 w Blade Sequence
Without breath-holding



T2 w Blade Sequence
Without swollen artifacts


Cao AAPM 2011

M Can we plan solely on MRI?

University of Michigan
Medical School

■ Issues:

- Bore size
- Distortion
- Electron density
- IGRT support



Cao AAPM 2011

M What sources of geometric errors and solutions are?

University of Michigan
Medical School

Sources of errors


- System level
 - B0 field inhomogeneity
 - Gradient non-linearity

Physics solutions

- System level
 - Better magnet design
 - On-line gradient distortion correction (GDC)
 - Algorithms to further correct any errors in system level

Cao AAPM 2011

Geometric Phantom to Map Homogeneity




Nina Hoven, Ulleval Hospital, Oslo, Norway

1.0T Philips panorama scanner

Distortion-free area:
Sagittal plane: 40 cm AP, 28 cm FH
Coronal plane: 34 cm FH, 36 cm LR
Transverse plane: 32 cm AP, 37 cm LR

Cao AAPM 2011

Gradient Distortion Correction




L Chen, Fox Chase Cancer Center, AAPM Summer School 2006
0.3 T scanner

Cao AAPM 2011

State-of-art 3T

Linear and high Orders correction




Cao AAPM 2011

Geometric distortion

| Sources of errors | Solutions |
|--|--|
| <ul style="list-style-type: none"> ■ Patient-level <ul style="list-style-type: none"> - Susceptibility - Fat/water chemical shift effect - Field strength - k-space trajectory - Gradient band width - Region: air, tissue, & bone interface | <ul style="list-style-type: none"> ■ Patient-level <ul style="list-style-type: none"> - Solutions: <ul style="list-style-type: none"> ■ B0 mapping ■ Rectification ■ published 15-20 years ago - Sub-mm for small FOV and 1-2 mm for large FOV distortions for SE and GE |

Cao AAPM 2011

Geometric Accuracy in Brain



Gradient Echo T1W Images from a 3T scanner


Registered to CT by rigid body transformation

Both superior and inferior portions of brain MRI are well registered to CT

Cao AAPM 2011

How can you get electron density from MRI?

- MR-CT alignment
- Atlas-based density insertion
- MR segmentation:
 - UTE imaging – attempts to directly visualize bone
 - Pattern learning to select candidate bone (versus air) features
- Hybrid approaches




G Bydder UCSD

Cao AAPM 2011

MRI-based patient modeling for RT planning

- Careful consideration of contrasts in MRI and human models permits image analysis to support:
 - Segmentation
 - Dose calculation
 - Image guided positioning



Cao AAPM 2011

Molecular/Functional/Metabolic MRI

- Molecular, metabolic and functional imaging
 - ^1H , ^{31}P and ^{13}C spectroscopy imaging
 - DCE and DSC imaging
 - DW and DT imaging
 - Other contrast agents, e.g., SPIO, Eovist
- Location, delineation, characterization, assessment of tumors and normal organs

Cao AAPM 2011




How to validate these imaging techniques for target definition?

- Pathological validation
 - Pathological specimen may not be easily obtained for certain organs, e.g., brain
- Pattern failure
 - Comparing the pattern pre RT with the recurrent pattern
- Prognostic and predictive factors
 - Via assessment of response or outcome to determine the subvolume of the tumor

Cao AAPM 2011



Primary Brain Tumor: GBM



CT

Post-Gd T1W

FLAIR


Underestimation

Overestimation/
Underestimation

Cao AAPM 2011



Proton Spectroscopy Imaging in Glioma


Metabolic Abnormality: CNI: Cho/NAA ≥ 2.0 SD

Cao AAPM 2011

Chan, J Neurosurg. 101:467,2004



Cho/NAA Abnormality in GBM







pre RT

4 months post RT

6 months post RT

LAPRIE, Int J Rad
Onc Biol Phys,
70:773-781, 2008

Cao AAPM 2011



Pathological Validation of 3D SI for Prostatic Tumor Localization

- UCSF study in 1999
 - 53 patients with biopsy-proved prostate cancer and subsequent radical prostatectomy with step-section histopathologic examination
 - T2W MRI:
 - sensitivity (77% and 81%), specificity (61% & 46%)
 - 3D MRSI (cho+Cr/citrate \geq 3SD):
 - sensitivity (63%) specificity (75%)
 - MRI+3D MRSI:
 - sensitivity (95% either test), specificity (91%)

Cao AAPM 2011

Validation of DCE and MRSI

- Schmuecking et al, Int J Radiat Bio 2009
- Evaluate quantitative DCE MRI and 1H MRS for the detection of prostate cancers and the delineation of intra-prostate sub-volumes for IMRT
- Groenendaal et al, Int J Rad Onc Biol Phys, 2010


Cao AAPM 2011

Delineation of Prostatic Cancers By DCE and MRS

- Schmuecking's study in 2009
 - Comparing quantitative DCE MRI and 1H MRS with these intraprostatic subvolumes with histology and cytokeratin-positive areas in prostatectomy species
 - DCE MRI: (1) 82% of sensitivity and 89% of specificity for localization of prostate cancers in left, right or both lobes; (2) able to detect the lesions > 3mm and/or containing >30% tumor cells; (3) similar to choline PET/CT
 - 1HMRS: (1) 55%-68% for sensitivity and 62%-67% for specificity; (2) able to detect the lesions > 8mm and/or containing >50% tumor cells


Cao AAPM 2011

DCE MRI Detection



Cao AAPM 2011

M DCE MRI vs MRS Detection



Cao AAPM 2011

M Challenges of MRI for RT

- Electron density
 - UTE MR imaging for bony structures
- Geometric accuracy
 - System level
 - patient-specific
 - Basic pulse sequences, e.g., GE and SE
 - QA/QC procedures
- Choose spatial resolution and plane orientation
- Position patients in the configuration of RT

Cao RSNA 2009

M Challenges of MRI for RT

- Sensitivity and specificity of each contrast or multiple contrasts for tumor delineation
- Reproducibility and uncertainty of metabolic and functional imaging
 - Spatial and amplitude
- Robustness of some of metabolic and functional imaging
- Optimize contrasts
 - Tumor specific
 - Optimal combinations of contrasts

Cao AAPM 2011



The Renaissance™
System 1000




Not Approved for Human or Animal Use

Cao RSNA 2009

University of Michigan Medical School

Localization of Intra-prostatic cancers by DWI



T2 WI Inversed DWI Pathologic Specimen

Positive detection rates of 6 observers:
 42–73% on T2WI alone
 58–80% on T2WI plus DWI


KAJIHARA, Int J Rad Onc Biol Phys, 74:399-403,2009

Cao RSNA 2010

University of Michigan Medical School

Delineation of Prostate GTVs

- Groenendaal's study
 - Comparing the GTVs delineated on DW and DCE MRI by a rad oncologist with the lesions (22) on prostatectomy specimens by a pathologist
 - 5 dominant intraprostatic lesions (>1cc) and 4 small lesions (>0.56 cc) detected by the Rad Oncologist based upon MRI
 - MRI GTVs of 5 DIL cover 44-76% of pathological tumor volumes but have 62-174% of the pathological tumor volumes



Cao RSNA 2010

University of Michigan Medical School


Delineation of Prostate GTVs

- Groenendaal's study
 - Sources of errors
 - Registration
 - Mis-matched characteristics between DW and DCE MRI (3 DIL), and negative on both DW and DCE MRI (1 DIL)
 - Solution
 - add 5 mm margin to the MRI-GTVs to improve the tumor volume coverage
 - The MRI-GTVs are 2.5-3 times as large as the pathological tumor volumes

Cao RSNA 2010

University of Michigan Medical School

DCE MRI– boost target?





Poorly perfused sub-tumor volume

Wang, Eisbruch, Cao, AAPM 2009

Cao RSNA 2010


M Gradient Distortion Correction

University of Michigan
Medical School



L Chen, Fox Chase Cancer Center,
AAPM Summer School 2006
0.3 T scanner


Cao AAPM 2011




State-of-art 3T scanner

M Soft Tissue Differentiation

University of Michigan
Medical School



CT




Post-Gd T1WI

Brain metastasis for SRS

Cao AAPM 2011

M CNI Abnormality vs Target in GKS of Recurrent GBM

University of Michigan
Medical School



>50%
Overlap
Survival:
15.7 m


<50%
Overlap
Survival:
10.4 m

Chan, J Neurosurg. 101:467,2004

M DCE MRI and MRS Detection

University of Michigan
Medical School

Schmuecking et al, Int J Radiat Bio 2009



The small lesion was missed by both DCE MRI and 1H MRS

Cao AAPM 2011