



## Personnel Protection During Fluoroscopic Procedures

Beth Schueler  
Mayo Clinic  
Rochester, Minnesota

## Medical Physicist's Role


- Procedure room design
- Movable and personal shielding device selection
- Fluoroscopic equipment selection
- Fluoroscopic system optimization
- Personnel exposure monitoring
- Personnel exposure investigation
- Radiation safety training

## Learning Objectives

- What is the nature of stray radiation in fluoroscopic procedures?
- What can be done to reduce occupational exposure?
- What regulations must be followed regarding personnel radiation exposure?

## Stray Radiation Sources

- X-ray tube leakage
- Scatter from patient and objects exposed to the primary beam




## X-Ray Tube Leakage

- Part 1020.30 (k) Leakage radiation from the diagnostic source assembly:
  - not to exceed 100 mR/hr (0.88 mGy/hr) at 1 m
  - measured at maximum kVp (125-150 kVp) and maximum continuous mA (3-5 mA)



## X-Ray Tube Leakage

- Leakage radiation is negligible compared to scatter
- typical amount of leakage radiation for fluoroscopy of an average adult abdomen: 5  $\mu$ Gy/hr at 1 m




## Scatter

- Number of scattered x-rays depends on
  - the patient entrance skin exposure rate
  - the patient entrance skin area
  - the beam energy
- Typical scatter value
  - 0.1% of incident exposure at 1 m for a 23-cm FOV
  - for an entrance skin exposure rate of 20 mGy/min, scatter at 1 m is 20  $\mu$ Gy/min or 1.2 mGy/hr

7

## Scatter vs Field Size


- For constant entrance exposure rate and constant kVp, scatter is proportional to the exposure field size



8

## Scatter vs Entrance Exposure


- For constant field size and constant kVp, scatter is proportional to the entrance skin exposure rate



9

## Scatter vs DAP


- For constant kVp, scatter is proportional to the dose-area-product rate (DAP)



10

## Scatter vs kVp


- For constant field size and constant entrance skin exposure rate, scatter increases with kVp



11

## Scatter: Distribution

- Scatter intensity is higher on the entrance side of the exposed volume
  - primary beam most intense on the entrance side
  - forward scatter is heavily attenuated




12

## Scatter: Distribution


- The entrance skin port should be considered the major radiation source for occupational exposure
- For actual imaging situations, the scatter distribution is also affected by
  - attenuation from unexposed tissue to the side
  - attenuation from the table and pad
  - shielding by image intensifier, equipment supports, lead shielding devices

13

## Stray Radiation: Horizontal Profile



## Stray Radiation: Vertical Profile



## Learning Objectives

- What is the nature of stray radiation in fluoroscopic procedures?
- What can be done to reduce occupational exposure?
- What regulations must be followed regarding personnel radiation exposure?

16

## Occupational Exposure Reduction

- What can operators do?
- What can assisting personnel do?
- What can medical physicists do?

17

## Radiation Safety: Operator

- Minimize fluoroscopy beam-on time
  - use short taps of fluoro instead of continuous exposure
  - use 5-min fluoro-on time warning to maintain awareness of total fluoroscopy time
  - utilize last-image-hold for image study, discussion and teaching
  - develop workflow procedures with ancillary personnel so that there is no fluoroscopic exposure when they need to be close to the patient

18

## Radiation Safety: Operator

- Increase distance from the patient
  - GI: “spoon” palpation device
  - vertebroplasty: PMMA injection device
  - other interventional devices: forceps



19

## Radiation Safety: Operator

- Remove the grid when possible
  - reduces patient entrance skin exposure and scatter by about a factor of 2
  - reduction in image contrast is minimal for
    - small patients and body parts
    - procedures requiring a large separation between the patient and image intensifier


20

## Radiation Safety: Operator


- Position yourself on the x-ray beam exit side of the patient when possible
  - Scatter intensity is lower on the image intensifier side as compared to the x-ray tube side
  - For C-arm fluoroscopy,
    - frontal projections: place the x-ray tube under the table
    - lateral and oblique projections: place the x-ray tube on opposite side of the patient

21

## Lateral Projection: Tube by Operator




## Lateral Projection: Tube Opposite Operator






24

## Radiation Safety: Operator

- Collimate tightly to the area of interest
  - reduces the entrance skin exposure area
  - improves image contrast
  - entrance skin exposure may increase
  - scatter will decrease with collimation as long as there is a decrease in DAP rate
  - similar result for wedge-shaped beam filters




- ### Radiation Safety: Operator
- Use optimal imaging chain geometry
    - General imaging guidelines
      - place the x-ray tube as far as possible from the patient
      - place the II as close as possible to the patient
      - use magnification modes sparingly
    - Adherence to these guidelines
      - reduces patient entrance skin exposure
      - increases spatial resolution
      - decreases scatter, but only when collimation is used
  - Use DAP rate as a guide to estimate scatter




- ### Radiation Safety: Operator
- Use movable shielding devices
    - ceiling-suspended shield
    - tableside shield
    - undertable shield
-

### Scatter Levels without Shielding



### Scatter Levels with Undertable Shielding



### Occupational Exposure Reduction

- What can operators do?
- What can assisting personnel do?
- What can medical physicists do?

33

### Radiation Safety: Assisting Personnel

- Make sure you are aware of beam-on cues
  - warning lights, audio warnings
- Increase your distance from the patient when possible
- Use movable shielding devices




34

### Radiation Safety: Operator and Assisting Personnel


- Be aware that scatter is significantly higher for large patients as compared to small patients
  - entrance skin exposure rate and kVp increases, resulting in higher scatter rates
  - image quality decreases, beam on-time may also increase

35


### Average-Size Patient



### Moderately Large-Size Patient



### Large-Size Patient




### Radiation Safety: Operator and Assisting Personnel

- Select appropriate personal shielding devices
  - Protective aprons
    - lead equivalent thickness of 0.5 mm recommended, required in some states
    - wrap-around style needed if side or back faces exposed patient volume
    - light-weight lead replacements (composites of materials such as barium, tungsten) reduce total weight by 20-30%

39

### Personal Shielding Devices

- Proper fit is critical
  - if 2 layers are needed to meet 0.5 mm lead-equivalent requirement, make sure overlap is sufficient
  - arm holes should not be too wide
  - neck line should not fall too low



40

### Personal Shielding Devices

- Thyroid shields
  - use recommended if monthly collar badge reading exceeds 4 mSv
  - significantly reduces effective dose with minimal inconvenience
- Leaded glasses
  - use recommended if monthly collar badge reading exceeds 4 mSv
  - side shields needed since user does not typically face the exposed patient volume

41

### Personal Shielding Devices

- Leaded gloves
  - 0.5 mm lead equivalent shielding recommended if hands must be in the primary beam
  - leaded surgical gloves provide partial shielding (30-40% attenuation) if hands are in high scatter area
    - provide minimal protection if hands are in the primary beam: ABC drives up dose rate, protection on one side only

42

## Occupational Exposure Reduction

- What can operators do?
- What can assisting personnel do?
- What can medical physicists do?

43

## Radiation Safety: Physicist

- When constructing or remodeling a fluoroscopic suite, ensure appropriate room design features are incorporated
  - ample procedure room size
  - ample control room size
  - large viewing windows
  - intercom system between procedure room and control room
  - patient table and control area positioned so that no straight-line scatter reaches

44

## Room Design Features

- Operator should be able to view all procedure room entrances from expected operating position
- Adequate structural shielding in walls, doors, and windows
  - consider structural shielding in procedure rooms where mobile C-arms are used routinely
- Movable shielding devices
  - ceiling-mounted shields should be included with fluoroscopic equipment installation


45

## Radiation Safety: Physicist


- Ensure fluoroscopic equipment radiation safety design features are considered when specifying a system for purchase
  - R&F systems: undertable x-ray tube configuration results in lower personnel exposure as compared to an overtable x-ray tube configuration (for non-remote applications)

46

### R&F System: Undertable X-Ray Tube



### R&F System: Overtable X-Ray Tube





## Equipment Design Features

- Undertable x-ray tube R&F systems (with lead II tower drapes attached) are preferable to C-arm configurations
- For biplane systems, a lateral C-arm that can be positioned with the II toward the operator is desirable



49

## Equipment Design Features

- Imaging chain components that allow positioning for optimal geometry
  - II can be moved close to the patient
  - eliminate spot-film device if not needed
  - x-ray tube located far from the patient
- Removable or retractable grid
- Last-image hold and fluoroscopy frame grab
- Graphic collimator and wedge filter position indicators on last-image hold

50

## Equipment Design Features

- Low dose fluoroscopy settings available
  - pulsed fluoroscopy at 15, 7.5 and lower frames per second
  - image processing: video frame averaging
  - use of spectral beam filtration
- Variable rate digital image acquisition available
- Patient dose indication available to the operator

51

## Equipment Design Features

- Digital acquisition zoom
- Low attenuation tabletop and pad
- Visible exposure warning lights
- Method to disable the fluoroscopy foot pedal to eliminate inadvertent exposure during patient removal or prep

52

## Radiation Safety: Physicist

- Implement and maintain a comprehensive quality control program
  - work with vendor representatives and operators to find optimal image processing and dose settings to meet clinical imaging needs
  - perform regular fluoroscopic equipment performance and radiation safety evaluation
  - perform personal protection device evaluation

53

## Learning Objectives

- What is the nature of stray radiation in fluoroscopic procedures?
- What can be done to reduce occupational exposure?
- What regulations must be followed regarding personnel radiation exposure?

54

## Maximum Permissible Doses

- Effective dose (whole body)
  - Annual: 50 mSv
  - Cumulative: 10 mSv × age in years
- Equivalent dose (partial body)
  - Lens of eye: 150 mSv
  - Skin, hands, feet: 500 mSv
- Embryo-fetus: 0.5 mSv/month

55

## Personal Dosimetry

- Single- and dual-monitor methods
  - single: above protective apron at collar level
  - dual:
    - Monitor 1. above protective apron at collar level
    - Monitor 2. under protective apron at waist level
- For the dual-monitor method, badges can easily be inadvertently switched
  - color coding of badges is helpful

56

## Dose Assignment

- Individual state regulations vary
- CRCPD Suggested State Regulations:
  - for single-monitor,  $H_E = \text{collar reading}$ , or
  - if single-monitor reading > 25% of the dose limit,  $H_E = 0.3 \times \text{collar reading}$ , or
  - for dual-monitor,  $H_E = 1.5 \times \text{waist reading} + 0.04 \times \text{collar reading}$ 
    - $H_E = \text{effective dose equivalent}$

57

## Dose Assignment

- NCRP Report No. 122 “Use of Personal Monitors to Estimate Effective Dose Equivalent and Effective Dose to Workers for External Exposure to Low-LET Radiation” 1995
  - for single-monitor,  $E = \text{deep dose reading}/21$
  - for dual-monitor,  $E = 0.5 \times \text{waist reading} + 0.025 \times \text{collar reading}$ 
    - $E = \text{effective dose}$

58

## Pregnant Workers

- Dose limits apply only after written declaration of pregnancy
- Conceptus dose assignment method
  - varies with individual state regulations
  - common method is to wear a monitor under the apron at the waist level
- Prediction of conceptus dose
  - estimate projected conceptus dose as 10% of collar badge for 0.5 mm lead-equivalent apron
    - Reference: Wagner and Hayman, Radiology 1982; 143: 559

59

## Pregnant Workers

- Prediction of conceptus dose
  - if pregnancy planned, monitor under-apron abdominal exposure prior to pregnancy
- Protection
  - 0.5 mm lead-equivalent aprons are sufficient
  - wrap-around styles critical if back may face radiation source
  - specially-designed maternity aprons available

60

## References

- American Association of Physicists in Medicine, "Managing the Use of Fluoroscopy in Medical Institutions." AAPM Report No. 58. Madison, WI: Medical Physics Publishing, 1998.
- Balter, S. Interventional Fluoroscopy. John Wiley and Sons (2001).
- Balter S. "Staff Safety in the Interventional Fluoroscopic Laboratory" in Intravascular Brachytherapy – Fluoroscopically Guided Interventions. Stephen Balter, Rosanna C. Chan, Thomas B. Shope, Jr (eds.). AAPM Monograph No. 28. Proceedings of the 2002 Summer School. Madison, WI: Medical Physics Publishing, pp. 137-163, 2002.
- Brateman L. "Radiation safety considerations for diagnostic radiology personnel." Radiographics 1999; 19:1037-1055.
- Geise RA and Hunter DW. "Personnel exposure during fluoroscopy procedures." Postgraduate Radiology 1988; 8:162-173.
- International Commission on Radiological Protection. Avoidance of Radiation Injuries from Medical Interventional Procedures. ICRP Publication 85. Ann ICRP 2000; 20(2).
- Marshall NW and Faulkner K. "The dependence of the scattered radiation dose to personnel on technique factors in diagnostic radiology." Med Phys 1996; 23:1271-1276.
- National Council on Radiological Protection and Measurements (NCRP). NCRP Report No. 105. "Radiation Protection for Medical and Allied Health Personnel." Bethesda, MD: NCRP, 1989.
- National Council on Radiological Protection and Measurements (NCRP). NCRP Report No. 122. "Use of Personal Monitors to Estimate Effective Dose Equivalent and Effective Dose to Workers for External Exposure to Low-LET Radiation." Bethesda, MD: NCRP, 1995.
- Wagner LK and Archer BA. Minimizing Risks from Fluoroscopic X-rays – Bioeffects, Instrumentation and Examination (2nd ed.). Houston, TX: Partners in Radiation Management (1998).
- Wagner LK and Hayman LA. "Pregnancy and women radiologists." Radiology 1982; 145:559-562.
- World Health Organization. Efficacy and Radiation Safety in Interventional Radiology. Geneva (2000).