

# Physics Uncertainty Issues for Prostate and Lung Cancers

Lei Dong, Ph.D.

Dept. of Radiation Physics


University of Texas M.D. Anderson Cancer  
Center, Houston, Texas

Symposium on the Promises and  
Perils of Proton Radiotherapy,  
Baltimore, MD, May 8-9, 2009

THE UNIVERSITY OF TEXAS  
MD ANDERSON  
CANCER CENTER  
*Making Cancer History*<sup>®</sup>

# Outline

- Review physics-related issues in proton therapy of prostate and lung cancers
- Discuss strategies and challenges for handling motion issue




# Potential Factors Affecting Prostate Treatment

- Target Motion (Prostate and SVs)
  - Inter-fractional rigid movement
    - Bony registration
    - Soft tissue target registration
  - Inter-fractional deformation
  - Intra-fractional motion
- Variations in Nearby Structures
  - Rectum filling
  - Bladder filling
  - Femoral heads

# Inter-fractional Variations

CT0



Planning contours mapped to 24 in-room CTs

# Inter-fractional Variations

Simulation CT is a snapshot of patient's anatomy!


Planning CT






Planning contours mapped to 24 in-room CTs

# Inter-fractional Variations

Simulation CT is a snapshot of patient's anatomy!




Planning contours mapped to 24 in-room CTs



Bi-lateral Proton


Trofimov et al. IJROBP 69(2): 444-453. (2007)

# Robustness of IMPT for Multi-beam IMPT?







# Femur Rotation



# Day-to-day variation in femur position





Planning CT




Anatomy over  
20 CT imaging  
sessions

# One Example




# Proton Distal Dose Variations Due to Femur Rotation

a)



Original Plan


b)



Daily CT

# Can CBCT/CT-on-Rails based IGRT help?


Direct alignment of soft tissue  
target (the prostate)



- Not entirely


- IGRT ensures the position of the prostate
- It will not guarantee normal tissue outside the target
  - Relative motion of the target relative to bony structures: the compensator can be off if bone is not aligned

a)




Original Plan

b)



Daily CT

Immobilization is important  
for proton therapy!



# Summary of Challenges in Prostate Proton Therapy

## 1. Uncertainties:

- CT numbers and stopping powers ratios,
- Inter-fractional changes of normal tissues in the beam path,
- Intra-fractional organ motion,


## 2. Strategies

- Optimization of beam angle selection
- Better immobilization techniques
  - Leg position
  - Rectal balloon




# Factors Affecting Proton Therapy of Lung Cancer


Free breathing treatment




# Proton vs. IMRT (motion effect to dose distribution)



# PTV Concept will not work



ICRU 62




Proton Therapy

# Effect of aperture expansion


RC smear: 0 mm  
AP margin: 0 mm  
Setup error: 0 mm  
Tumor motion: 0 mm

PLANNED



RC smear: 0 mm  
AP margin: 0 mm  
Setup error: 10 mm  
Tumor motion: 0 mm

“DELIVERED”




Isodose  
levels  
20  
50  
80  
95  
100

# Effect of aperture expansion


RC smear: 0 mm  
AP margin: 5 mm  
Setup error: 0 mm  
Tumor motion: 0 mm

PLANNED



RC smear: 0 mm  
AP margin: 5 mm  
Setup error: 10 mm  
Tumor motion: 0 mm

“DELIVERED”




Isodose  
levels  
20  
50  
80  
95  
100

# Effect of aperture expansion


RC smear: 0 mm  
AP margin: 10 mm  
Setup error: 0 mm  
Tumor motion: 0 mm

RC smear: 0 mm  
AP margin: 10 mm  
Setup error: 10 mm  
Tumor motion: 0 mm

PLANNED



“DELIVERED”




Isodose  
levels  
20  
50  
80  
95  
100

Aperture expansion alone is not a solution

# Effect of range compensator smearing to ensure target coverage

- The bolus thickness at point P, as calculated by the 'simple' technique, is replaced by the thinnest bolus thickness calculated anywhere within  $\pm d$  of the point P




# Effect of range compensator smearing


RC smear: 5 mm  
AP margin: 0 mm  
Setup error: 0 mm  
Tumor motion: 0 mm

RC smear: 5 mm  
AP margin: 0 mm  
Setup error: 10 mm  
Tumor motion: 0 mm

PLANNED



“DELIVERED”



Isodose  
levels  
20  
50  
80  
95  
100




# Effect of range compensator smearing


RC smear: 10 mm  
AP margin: 0 mm  
Setup error: 0 mm  
Tumor motion: 0 mm

RC smear: 10 mm  
AP margin: 0 mm  
Setup error: 10 mm  
Tumor motion: 0 mm

PLANNED



“DELIVERED”




Isodose levels  
20  
50  
80  
95  
100

Smearing alone is not a solution

# Compensator Smearing and Aperture Expansion


RC smear: 0 mm  
AP margin: 0 mm  
Setup error: 0 mm  
Tumor motion: 0 mm

PLANNED



RC smear: 7.5 mm  
AP margin: 7.5 mm  
Setup error: 0 mm  
Tumor motion: 0 mm

PLANNED




Isodose  
levels  
20  
50  
80  
95  
100


# Compensator Smearing and Aperture Expansion

RC smear: 7.5 mm  
AP margin: 7.5 mm  
Setup error: 0 mm  
Tumor motion: 0 mm

PLANNED






RC smear: 7.5 mm  
AP margin: 7.5 mm  
Setup error: 5 mm




Isodose  
levels  
20  
50  
80  
95  
100

How to evaluate a plan in the presence of anticipated motion?




# CT Number uncertainties have a significant impact on proton dose distributions

But commonly it's not visible on proton plans




# Scatter plot of various human biological tissues listed in ICRU Report 44 and ICRP Report 23

Human biological tissues have different electron density ratios (HUs) and effective atomic numbers





Proton Stopping-Power-Ratio depends on the effective atomic numbers (~ 3 - 4 % uncertainties)



# Improving CT number Accuracy

- High Energy CT imaging to reduce metal artifacts
- Dual-energy CT imaging to characterize material compositions
  - Electron density
  - Effective atomic number



# Challenges in planning moving target

- 4D dose calculation
  - Accumulating dose distributions voxel-by-voxel and phase-by-phase
    - Require deformable image registration and re-computing plans for each breathing phase using 4D CT
- 4D treatment plan optimization
  - Automatically adjust distal margin, aperture, compensator smearing parameters to satisfy target coverage and normal tissue sparing
- 4D treatment plan evaluation

None of these are currently available in a commercial product!



# Proton Planning Guide for Non-Gated Treatment

## - A practical approach by MDACC and MGH

- Target delineation
  - ITVs delineated on a 4D CT set
- Dose calculation
  - Using the time-averaged CT based on the 4D CT set
  - Override iGTV with tissue density
    - A compromise to cover the moving target without altering TPS internal compensator design
- Planning
  - Evaluate DVHs using contours derived from Mid-ventilation CT
- Final Plan Summary
  - Copy the plan to Mid-ventilation CT without any density override
  - Report DVHs using lung volumes defined either for Mid-Ventilation or the Expiration phase

# Special Situations

- If gated treatment is not used

- Tumor near the diaphragm


- The patient is not a good candidate for proton therapy if any beam can not avoid a portion of the diaphragm before entering the treatment target.

- We don't have a good solution!


- Stomach filling

- If we cannot avoid the beam direction through the stomach, the patient is not a good case for proton therapy.

# Esophagus cancer near GE junction




Green: density changes



Red: GTV


# More problems: Distal Edge Degradation Caused By Heterogeneities

Heterogeneity




# “Texturized” Lung Equivalent Material Monte Carlo vs. TPS Dose Calculations

## Simulations with Validated Monte Carlo



- Texturized phantom:  $3 \times 3 \times 3 \text{ mm}^3$  voxels
- Random 20% with unit density, rest air
- Average density = 0.2 gm/cc

# Range Degradation Through Lung Equivalent Material - Homogeneous vs. "Texturized"



Calculation Error!

# Summary of Physics Uncertainty Issues

- Proton beams are more sensitive to
  - Organ motion
  - Anatomy changes in the beam path
  - CT number accuracy
  - Treatment devices in the beam path (for example, the couch, immobilization devices)
- IGRT and motion management may help
- Proton plans are difficult to evaluate
- Many challenges remain due to practical reasons
  - Compromises may have to be made for routine practice
    - **Use common sense!**