

Acceptance Testing and Commissioning of Monte Carlo Dose Calculation Systems

Bruce Curran
University of Michigan Medical Center
Ann Arbor, MI

Disclosure

Some of the work discussed in this talk was supported by NOMOS Corporation

Acknowledgements

- UCSF
 - Cynthia Chuang, Bruce Faddegon, Lynn Verhey
- McGill University
 - Emily Heath, Jan Seuntjens
- Chiba Medical Center
 - Yuichiro Narita
- Mayo Scottsdale
 - Gary Ezzell

Acceptance Testing

- Performing an agreed set of tests to establish that the delivered product meets specifications
 - Tests often pre-packaged by manufacturer
 - Scope of tests is often limited
 - Does the calculation run?
 - Often, there is no statement of accuracy in dose algorithms
 - Most of the work is in commissioning

Acceptance Testing

- What Acceptance Testing does do
 - Verifies that hardware / software is installed correctly
 - Gives user some training in operation of software
- What Acceptance Testing does NOT do
 - Verify accuracy of algorithm over a range of clinical situations
 - Verify accuracy of your specific beam data entry
 - Provide complete training on operation and evaluation of the dose calculation

Special Issues for Acceptance Testing for Monte Carlo Dose Calculation Systems

- In general, treatment planning systems only specify accuracy in implementing an algorithm
 - No direct ties to measurements
 - No requirement to compare with actual data during acceptance testing
- Monte Carlo systems may specify dose accuracy, not just algorithmic accuracy

Sample Criteria for Dose Algorithm Accuracy


Situation	Absolute Dose (%) ^a	Central Ray (%)	Inner Beam (%)	Penumbra (mm)	Outer Beam (%)	Build-up Region (%)
A. Homogeneous Phantoms						
Square fields	0.5	1	1.5	2	2	20
Rectangular fields	0.5	1.5	2	2	2	20
Asymmetric fields	1	2	3	2	3	20
Blocked fields	1	2	3	2	5	50
MLC-shaped fields	1	2	3	3	5	20
Wedge fields	2	2	5	3	5	50
External surface variations	0.5	1	3	2	5	20
SSD variations	1	1	1.5	2	2	40
B. Inhomogeneous Phantoms^{b,c}						
Slab inhomogeneities	3	3	5	5	5	-
3-D inhomogeneities	5	5	7	7	7	-

From Van Dyk Presentation (AAPM-2004) [adopted from IAEA TRS430 and AAPM TG-53]

Special Issues for Acceptance Testing for Monte Carlo Dose Calculation Systems

- Head Models
 - Often compartmentalized (e.g. hidden) from user
 - Do they have the correct information for your machine?
 - Proprietary information
 - Often simplified
 - What details are necessary to accurately model your treatment delivery system?
 - May be beyond the scope of a user to define
 - Complex geometries / materials: Should each user be required to define the same class of delivery device?

DYNVMLC Component Module



As a result acceptance testing and commissioning cannot truly be separated for Monte Carlo Dose Calculation Systems

What can we expect today?

- References to Geometric Models
 - May include proprietary documents
- Sample data showing clinically relevant results
 - Same beam energy, ...
 - Data should include heterogeneous materials
 - Most useful if comparable measurements shown for accuracy

Monte Carlo Systems: Raising the Bar?

Most medical physicists expect (?require?) Monte Carlo based dose calculations to be

- more accurate
- more precise
- more reliable

Measurements for Monte Carlo verification are

- affected more strongly by detector characteristics
- significantly subject to variations in protocol (TG-21, TG-51)

Monte Carlo Dose Accuracy?

- Dose Accuracy
 - 2% under all conditions?
 - (or just high dose / low gradient areas?)
 - 2 mm in high gradient regions
 - Is that good enough in build-up regions?
 - ?? in low dose / low gradient regions
 - under MLCs, jaws, blocks

Issues in verifying MC Dose Accuracy

- CT Issues
 - Mass Density, Material Calibration
 - Electron Density Errors
 - Edge Effects
- Dose Specification
 - to water or to media?
- Accelerator Characteristics
 - Class variations
- Measurement Equipment
- Training

Edge Effects (Ideal Phantom)


-1000	-1000	-1000	-1000	-1000	-1000	-1000	-1000	-1000	-1000
-1000	-1000	-1000	-1000	-1000	-1000	-1000	-1000	-1000	-1000
-1000	-1000	-1000	-1000	-1000	-1000	-1000	-1000	-1000	-1000
0	0	0	0	0	0	0	0	0	0
0	0	0	0	0	0	0	0	0	0
0	0	0	0	0	0	0	0	0	0
0	0	0	0	0	0	0	0	0	0
0	0	0	0	0	0	0	0	0	0
0	0	0	0	0	0	0	0	0	0
0	0	0	0	0	0	0	0	0	0
0	0	0	0	0	0	0	0	0	0

Edge Effects (Scanned Phantom)

-1000	-1000	-1000	-1000	-1000	-1000	-1000	-1000	-1000	-1000
-400	-400	-400	-400	-400	-400	-400	-400	-400	-400
-200	-200	-200	-200	-200	-200	-200	-200	-200	-200
0	0	0	0	0	0	0	0	0	0
0	0	0	0	0	0	0	0	0	0
0	0	0	0	0	0	0	0	0	0
0	0	0	0	0	0	0	0	0	0
0	0	0	0	0	0	0	0	0	0
0	0	0	0	0	0	0	0	0	0
0	0	0	0	0	0	0	0	0	0
0	0	0	0	0	0	0	0	0	0


Edge Effects

Imported Solid Water CT Scan



Edge Effects

Imported Solid Water CT Scan



Edge Effects

- Caused by sudden, sharp edges in scanned volume
- Not just interpolation effect
- Can result in 'extending' the phantom several mm beyond its actual boundary
- Can offset depth dose, build-up regions

Electron Density

(Actual Phantom)

0	0	0	0	0	0	0	0	0	0
0	0	0	0	0	0	0	0	0	0
0	0	0	0	0	0	0	0	0	0
1150	1150	1150	1150	1150	1150	1150	1150	1150	1150
1150	1150	1150	1150	1150	1150	1150	1150	1150	1150
1150	1150	1150	1150	1150	1150	1150	1150	1150	1150
0	0	0	0	0	0	0	0	0	0
0	0	0	0	0	0	0	0	0	0
0	0	0	0	0	0	0	0	0	0
0	0	0	0	0	0	0	0	0	0

Electron Density (Scanned Phantom)

0	0	0	0	0	0	0	0	0	0
0	0	0	0	0	0	0	0	0	0
0	0	0	0	0	0	0	0	0	0
1150	1150	1000	900	650	650	900	1000	1150	1150
1150	1150	1000	900	650	650	900	1000	1150	1150
1150	1150	1000	900	650	650	900	1000	1150	1150
0	0	0	0	0	0	0	0	0	0
0	0	0	0	0	0	0	0	0	0
0	0	0	0	0	0	0	0	0	0
0	0	0	0	0	0	0	0	0	0

Electron Density

- Caused by Beam Hardening Correction of CT Scanner
- ‘Corrects’ for higher attenuation in center of patient
- Primarily affects non-patient geometries, e.g. rectangular solids, slab heterogeneities

Dose Specification


- Most centers calibrate dose to water
- MC calculations often specify dose according to media of voxel
 - CT to Material Calibration Curve
- Can result in differences of 1-3% depending upon region being measured

Accelerator Characteristics


(Lots of Questions)

- What physical parameters vary from machine to machine within an accelerator class?
 - Energy
 - Spectrum
 - Beam spot size
 - Energy / Intensity Distribution
- How can we handle manufacturing variations?
- How does a manufacturer compensate for measurement variations?


Effect of electron energy




Energy Determination




Energy Extraction




Measurement Variations




Off-axis Profiles (d=1.5 cm)




Ionization Chamber Volume Averaging



Low-Dose Results




Buildup region (10x10 cm²)



Training

- Proper Measurement Techniques
 - Ion Chamber
 - Film
- Isodose Interpretation
 - Statistical Variation vs Imprecision

What are we drawing?



Summary

- Acceptance Testing of Monte Carlo Dose Calculation Systems can require significantly more effort than traditional dose algorithms
- Imperfections in both CT and Measuring Systems can result in imperfect data for use in Monte Carlo System commissioning
- MC Dose Calculation System Commissioning will require efforts by both manufacturer and user
- There will need to be some re-thinking of our traditional planning processes in order to fully utilize Monte Carlo techniques.