

VCU *The Department of Radiation Oncology*

Monte Carlo II Planning:

Applications in IMRT QA, IMRT Optimization, Motion Compensation, and 4D Dose Calculation

J.V. Siebers
Virginia Commonwealth University
Medical College of Virginia Hospitals
Richmond, Virginia USA

2006 AAPM Continuing Education

Objectives

- To understand methods and role of MC for
 - IMRT dose calculation
 - Patient-specific IMRT QA
 - IMRT optimization
 - 4D dose calculation

2006 AAPM MC Continuing Education

What is new

- Vendors are showing IMRT specific IMRT products
 - CMS : Monaco
 - Elekta
- MC optimized plans are being achieved in <30 minutes
- MC is being applied to fluence prediction with conventional algorithms for patient dose calculation

2006 AAPM MC Continuing Education

Justification of Monte Carlo for IMRT

- IMRT is so complex, patient specific QA is mandated
 - Planning dose prediction
 - Information transfer
 - Device delivery
- Conventional dose algorithms can be inaccurate for
 - Small fields
 - Regions of dose gradients (radiation disequilibrium)
 - Heterogeneous conditions
- If algorithm was always accurate, there would be no need for per-patient algorithm QA

2006 AAPM MC Continuing Education

Dose Calculation Basics


2006 AAPM MC Continuing Education


Dose Calculation is a Two Step Process

- Incident fluence prediction
- Energy deposition in the patient/phantom

2006 AAPM MC Continuing Education

MC can be used for both steps, or only one step of the process






2006 AAPM MC Continuing Education

Intensity Modulation Matrix for IMRT Fluence Prediction


Conventional algorithm use of intensity matrix


$$\Psi(x, y)_f = \Psi(x, y)_i \times I(x, y)$$



MC utilizing TPS intensity matrix


$$w_f(x, y) = w_i(x, y) \times I(x, y)$$






2006 AAPM MC Continuing Education


Caution: MC using TPS Intensity Matrix does NOT improve fluence prediction




Convolution



Monte Carlo (same fluence)




Analytic MLC intensity algorithms approximate effects of MLC scatter, beam hardening, ...




2006 AAPM MC Continuing Education

MC for Fluence Prediction

Direct particle transport



- Simulate individual particles through MLC
- Geometric details can be fully included
- Physics (leakage, scatter) inherently accounted
- SMLC and DMLC





2006 AAPM MC Continuing Education

MC for Fluence Prediction

Varian 120 leaf

- General purpose MC algorithms (MCNP, GEANT, ...)
- Accelerator simulation algorithms (BEAM, BEAMnrc)
- Specialized algorithms
 - Fast ray tracing
 - Efficient particle use
 - Simplified geometry
 - Simplified physics






using Education

Efficient particle usage


Average attenuation from multiple time samples



Determine mean weight by sampling at multiple random "times".

$$w_f = \frac{w_i}{N} \sum_{k=1}^N e^{-\mu(E)t_k / \cos \theta_z}$$



Methods for MC For Patient Dose Prediction



2006 AAPM MC Continuing Education

MC for patient dose prediction Bixel based MC calculations



- Incident fluence sub-divided into ~1x1 cm² beamlets
- One-time dose computation for each beamlet to generate bixels
- Optimization determines bixel weights

2006 AAPM MC Continuing Education

MC for patient dose prediction Bixel based MC calculations



- Advantages:
 - One-time dose computation
- Disadvantages:
 - Limited resolution
 - Memory storage
 - Neglects effect of beam delivery

2006 AAPM MC Continuing Education

MC for patient dose prediction Field-based dose calculation



- Compute dose for entire beam
- For SMLC, can compute per-segment dose computation (use in segment weight optimization)
- Advantages:
 - Fluence can include MLC effects (leakage / scatter)
- Disadvantage:
 - Re-computation each time fluence changes (different leaf positions) optimization)

2006 AAPM MC Continuing

MC for patient dose prediction Field-based dose calculation

- Compute dose for entire beam
- For SMLC, can compute per-segment dose computation (use in segment weight optimization)
- Advantages:
 - Fluence can include MLC effects (leakage / scatter)
- Disadvantage:
 - Re-computation each time fluence changes (different leaf positions) optimization)








2006 AAPM MC Continuing

MC for patient specific IMRT QA




2006 AAPM MC Continuing Education



MC during IMRT optimization Pre-optimization

- Optimization will converge in fewer iterations if a good initial guess is provided to the optimizer
- MC optimization should be preceded with pre-optimization using faster/approximate algorithms

2006 AAPM MC Continuing Education



MC IMRT optimization (a la Hyperion)

- PB pre-optimization
- MC computed leaf sequences (step-and-shoot segments)
- Optimize segment weights
 - Account for heterogeneities
 - Accounts for leaf effects
 - Adjust segments based on PB gradients

2006 AAPM MC Continuing Education

DMLC Field-based optimization

Requires including MLC leaf sequences into MC optimization

```

    graph TD
      1[Initial Intensity I(x,y)] --> 7[Create Leaf Sequence]
      7 --> 8[Create Deliverable Intensities I_d(x,y)]
      8 --> 2[Compute Dose D_c]
      2 --> 3[Evaluate Plan Objective]
      3 --> 4{Converged?}
      4 -- No --> 7
      4 -- Yes --> 5[Optimized Intensity I_o(x,y) and Dose D_o = I_o]
      5 --> 6[Output I_o(x,y)]
      2 -.-> Note1[MC_MLC - new leaf sequence at each iteration]
      5 -.-> Note2[Final dose is deliverable]
    
```

2006 AAPM MC Continuing Education

Clinical Case Comparisons

2006 AAPM MC Continuing Education

Patient Comparison

Analytic-Analytic SC optimized vs MC^{AnalyticLeafSequences}-MC_{MLC} Recomputed with MC_{MLC}

66 Gy 60 Gy 57 Gy 40 Gy 20 Gy

66 Gy Hot-Spot location

57 Gy line not cover PTV

2006 AAPM MC Continuing Education

Optimized with MC

(a) Approved plan that did not agree with MC

(b) MC optimized plan restores target coverage

Initial desired dose distribution was achievable, but it required different intensities / leaf sequences than predicted by SC to be achieved in the patient

2006 AAPM MC Continuing Education




MC to reduce errors

Head and Neck IMRT plan

MC deliverable optimization can restore original optimized plan

DPE = Dose Prediction Error OCE = Optimization Convergence Error


2006 AAPM MC Continuing Education



MC optimization times

- Hyperion
 - Step-and-shoot
 - XVMC
 - < 1 hour
- VCU
 - Dynamic MLC field-based dose computation
 - VMC++ with hybrid algorithm
 - 20 processors
 - < 30 minutes (2 MC iterations)
 - Poster TU-EE-A2-04

2006 AAPM MC Continuing Education



Summary MC applications in IMRT


- MC can account for heterogeneities and/or fluence prediction errors
- Useful for patient specific QA
- MC can improve IMRT dose accuracy (QA implications)
- MC optimization has become clinically practical

Vendors demonstrating MC-IMRT products

2006 AAPM MC Continuing Education

Four-dimensional Monte Carlo dose calculation


- Account for intra- and inter-tissue motions in dose evaluation
- Most mobile tissues (lung) also have greatest benefit from MC



2006 AAPM MC Continuing Education

4D Dose Calculation Method

- Create deformation fields
- Define anatomy on reference CT
- Map anatomy to all CT sets
- Create treatment plan on reference CT
- Create/Apply treatment plan to other CT sets, Compute Dose
- Combine dose distributions and display on reference CT



2006 AAPM MC

MC 4D dose calc

Patient: Transport techniques

Either


- Per-anatomy (phase) rectilinear dose grid

or

- Deformed dose grid (Heath *et al*, Med Phys, 36 6)


Patient: Statistics

- Can compute to relatively poor statistics on each anatomic sample




2006 AAPM MC Continuing Education

Isodoses from 6% uncertainty calculation per beam at peak inhale



2006 AAPM MC Continuing Education


Map per-phase 6% dose distributions to reference (equi-time)




2006 AAPM MC Continuing Education

6% uncertainty beam

Combined 8 phases



2006 AAPM MC Continuing Education



4D MC Summary

- Advantages of MC for 4D dose evaluation
 - (1) higher accuracy for calculation in electronic disequilibrium conditions such as those encountered during lung radiotherapy
 - (2) if deformable image registration is used, the calculation time for Monte Carlo is **≈independent** of the number of 3D CT image sets constituting a 4D CT (same time for 4D and 3D calculation)

2006 AAPM MC Continuing Education

Thank You for Your Attention

- Note
 - Several groups are using MC and IMRT
 - 14 abstracts at this meeting on MC and IMRT
- Acknowledgements/Contributors
 - Paul Keall
 - Ivaylo Mihaylov
 - Weidong Li
 - Iwan Kawrakow
 - Radhe Mohan
 - S. Joshi

2006 AAPM MC Continuing Education