





Monte Carlo Dose Calculation for Clinical IMRT Treatment Planning

J.V. Siebers
*Virginia Commonwealth University
Medical College of Virginia Hospitals
Richmond, Virginia USA*
2008 AAPM Continuing Education

Funded by: NIH R01 CA98524



Objectives





- To understand stages of an MC-IMRT dose calculation
- To illustrate the role of MC for
 - Treatment Planning System QA
 - Patient-specific IMRT QA
 - Clinical IMRT plan optimization
- To compare clinical IMRT dose calculations/optimization with MC-based results


Virginia Commonwealth University

Effect of Heterogeneities on Dose

Errors in other algorithms are largest for small fields in heterogeneous geometries

Anfield et al., Med. Phys., 27 (6) 2000




Virginia



Justification of Monte Carlo for IMRT

- IMRT mandates patient specific QA
 - Treatment planning system dose prediction
 - Information transfer
 - Device delivery
- Conventional TPS dose algorithms can be inaccurate for
 - Small fields
 - Regions of dose gradients (radiation disequilibrium)
 - Heterogeneous conditions
- If dose prediction was always accurate, there would be no need for per-patient QA
 - Simpler methods exist for checking information transfer and device delivery than doing measurement-based QA


Virginia Commonwealth University



Particle Source

- Direct Simulation through tx head
 - On-the-fly simulation
 - Stored particles (phase space file)
- Source Models
 - Simulation-based models
 - Measurement-driven empirical source models

Covered in detail elsewhere (TG-105, ...)
 Overall accuracy is no better than least accurate sub-component

VCU Virginia Commonwealth University

Field-specific beam modifiers

- May be part of the source model

Or

- Explicit simulation

MLC Movie courtesy
 Jai-Woong Yoon
 Stanford


VCU Virginia Commonwealth University

Preliminary Requirements


- Integrated MC system
- Commissioned MC algorithm
 - Source model tuning
 - Open field matching
 - MLC test fields
 - IMRT test fields
- Other
 - See TG-105 ++

VCU Virginia Commonwealth University


MC for patient specific IMRT QA



VCU Virginia Commonwealth University



MC QA Results




	DMLC-VER	MCNUCHECK
Brainstem D.O. 01	4742.8	4848.4 (-2.23)
Brainstem D.O. 05	4604.4	4596.8 (+0.16)
Brainstem D.O. 50	3946.9	3936.8 (+0.09)
Brainstem D.O. 90	490.2	236.0 (+51.86)
Brainstem D.O. 98	432.3	181.8 (+57.96)
Cord D.O. 01	4049.2	4135.2 (-8.59)
Cord D.O. 05	3980.4	4104.1 (-3.11)
Cord D.O. 50	3713.8	3639.5 (+2.00)
Cord D.O. 90	3422.9	89.3 (+497.39)
Cord D.O. 98	1713.5	49.5 (+497.11)
Larynx D.O. 01	5180.5	5435.1 (-4.91)
Larynx D.O. 05	5072.3	5252.1 (-5.52)
Larynx D.O. 50	4634.8	4772.7 (-2.98)
Larynx D.O. 90	4222.1	4389.8 (-3.70)
Larynx D.O. 98	3988.1	4164.8 (-4.43)
LTParotid D.O. 01	5809.3	6036.9 (-3.94)
LTParotid D.O. 05	5634.9	5973.8 (-4.24)
LTParotid D.O. 50	1895.9	2180.9 (-15.03)
LTParotid D.O. 90	1260.7	1863.7 (-32.99)
LTParotid D.O. 98	995.8	1519.0 (-51.54)
Mandible D.O. 01	6182.7	6301.1 (-1.92)
Mandible D.O. 05	5584.7	5750.8 (-1.97)
Mandible D.O. 50	3524.3	3697.2 (-4.90)
Mandible D.O. 90	2575.0	2759.8 (-7.18)
Mandible D.O. 98	2168.3	2378.1 (-8.68)
PTV1 D.O. 01	7465.3	7607.0 (-1.90)
PTV1 D.O. 05	7422.9	7547.7 (-1.68)
PTV1 D.O. 50	7246.1	7336.9 (-1.25)
PTV1 D.O. 90	7093.6	7140.7 (-0.66)
PTV1 D.O. 98	6993.7	7036.8 (-0.65)
PTV2 D.O. 01	7437.9	7565.4 (-1.71)
PTV2 D.O. 05	7368.4	7482.9 (-1.55)
PTV2 D.O. 50	6899.4	6992.4 (-1.50)
PTV2 D.O. 90	6289.0	6338.6 (-0.79)
PTV2 D.O. 98	6092.8	6123.0 (-0.50)
PTV3 D.O. 01	7399.6	7518.5 (-1.61)
PTV3 D.O. 05	7281.7	7383.1 (-1.38)
PTV3 D.O. 50	5871.0	5956.2 (-1.45)
PTV3 D.O. 90	5567.1	5587.9 (-0.37)
PTV3 D.O. 98	5446.2	5483.6 (-0.69)
RightParotid D.O. 01	7162.6	7373.0 (-2.94)
RightParotid D.O. 05	6946.3	7143.7 (-1.80)
RightParotid D.O. 50	5442.6	5693.7 (-4.61)
RightParotid D.O. 90	3256.6	3465.9 (-6.23)
RightParotid D.O. 98	2776.0	2851.0 (-2.70)

VCU Virginia Commonwealth University


QA Summary

- Pinnacle vs MC Dose
- Pinnacle & VCU IMRT
- Any IMRT Tx Site
- $\Delta = -0.77\%$
- $\sigma = 1.43\%$
- 7% of cases have > 3% dose difference



VCU Virginia Commonwealth University


Head and Neck Example



VCU Virginia Commonwealth University

Head and Neck Example


- PTV $D_{98} \Delta = -2.8\%$
- LtParotid $D_{50} \Delta = -11\%$
- Cord $D_1 \Delta < 1\%$






VCU Virginia Commonwealth University



Head and Neck Example - Pinnacle using MC Fluence-

- PTV $D_{98} \Delta = -2\%$
- LtParotid $D_{50} \Delta = -3.5\%$



VCU Virginia Commonwealth University






MC during IMRT optimization


Pre-optimization

- **Since**
 Optimization will converge in fewer iterations if a good initial guess is provided to the optimizer
- **Then**
 MC optimization should be preceded with pre-optimization using faster/approximate algorithms


VCU Virginia Commonwealth University




Optimization Sequence




Optimization Sequence







Optimization Sequence




Optimization Sequence





Head and Neck



1 MC-based correction sufficient

MC vs One Correction Hybrid Comparisons

Structure	Dose Volume Constraint (Gy)	DPE of initial calc D^{95} (Gy)	DPE of D^{95} of PB (Gy)	OCE of initial optimization D^{95} (Gy)	OCE of D^{95} of PB (Gy)
GTV	$D_{95} > 68.1, 70.8, \text{ or } 73.8$		-0.05 ± 0.09		0.12 ± 0.17
CTV	$D_{95} > 60.0$		-0.06 ± 0.06		0.07 ± 0.12
ETV	$D_{95} < 54.0$		-0.02 ± 0.05		0.05 ± 0.09
Parotids	$D_{12} < 30.0$		-0.06 ± 0.04		0.00 ± 0.11
Cord	$D_{12} < 45.0$		-0.03 ± 0.07		-0.07 ± 0.18
Brainstem	$D_{12} < 55.0$		-0.00 ± 0.19		0.01 ± 0.06

Time Reduction Maximum tolerated σ_{Dose} during optimization

- 20 patients (10 Prostate, 10 Head & Neck)
- Optimize with different σ_{Dose}
 - MC optimization & One-MC correction optimization
 - $\sigma_{\text{Dose}} = 2, 5, 8, 10, 12, 15, 20$ % per beam
 - Achieved stats 0.9%-9% wrt local dose
- Verification calculation with 2% per beam statistics
- Compare with 2% per beam MC optimization
 - DVH indices used to compute PQS during optimization with statistical equivalence tests
 - Isodoses, Full DVH's,
- Rejection Criteria: 1 Gy Dose Difference ($\sim 1.5\%$ of Rx dose) between optimized plans

Statistical Equivalence Test Results

Prostate

- Dose window size (in cGy) required to state optimization results are equivalent to the MC 2%/beam optimization result

	2%*	5%*	8%*	10%*	12%*	15%*	20%*
MC	--	77	47	88	131	214	180
OC	40	63	51	67	130	195	181

Maximum tolerated < 100 cGy threshold.

10% per beam corresponds with per-plan $\sigma \sim 5\%$ wrt D_{local}

*Nominal per-beam statistical precision

Dose computation time


Statistical Precision (%)		Dose Calculation Time Required for Optimization (min)			
Nominal Per Beam	$\overline{\sigma}_{D_{max}}$ Achieved	$\overline{\sigma}_{D_{min}}$ Achieved	MC _{opt}	Hybrid	One Correction
2	0.9	0.6	700	326	164
5	2.4	1.5	102	48	25
8	3.7	2.4	40	18.6	9.6
10	4.7	3	25.7	12.1	6.2
12	5.6	3.5	17.6	8.3	4.3

- Scaled to single OLD 3 Ghz Xeon processor (Current CPUs 3X faster)
- 32-bit version of VMC++, DMLC, BEAM (64-bit version reported to be faster)
- Phase-space style source (source model would be 2x faster)
- 4x4x4 mm³ voxels


- Head/Neck optimizations take ~1.5x longer (split beams)

VCU Virginia Commonwealth University


MC Lung Case: QA



Lungs




Monte Carlo Optimization



- Optimize with MC
- MU Check (QA) with SC

VCU Virginia Commonwealth University



MC times

(OLD CPUs @ VCU, new CPUs 3X faster)





- **QA**
 - VCU
 - ✓ dosxyznrc: 4.1 hours (20*3 Ghz CPU, $\sigma_{Dmax}=0.6\%$)
 - ✓ dosxyznrc: 38 min (20*3 Ghz, $\sigma_{Dmax}=1.5\%$)
 - ✓ VMC++: 25 min (1*3 Ghz CPU, $\sigma_{Dmax}=1.5\%$)
 - Yamamoto (PMB 52 2007)
 - ✓ 15 hrs on 28 CPU cluster
- **Optimization**
 - Hyperion / CMS (Step-and-shoot)
 - ✓ XVMC, < 1 hour (multi-core, multi CPU system)
 - VCU (dynamic MLC)
 - ✓ Full MCopt with dosxyznrc: 12.5 hours (1 CPU, $\sigma_{Dmax}=3\%$)
 - ✓ One Correction Hybrid with VMC++: 6 min (Opt, $\sigma_{Dmax}=3\%$) + 25 min (final calc, $\sigma_{Dmax}=1.5\%$) (1 CPU)


VCU Virginia Commonwealth University

Summary

MC applications in IMRT

- MC can account for heterogeneities and/or fluence prediction
- MC useful for 3D patient specific QA
- MC optimization is clinically practical
 - Calculation times few minutes
 - Reduces normal structure dose



VCU Virginia Commonwealth University

Thank You for Your Attention

- **Acknowledgements / Contributors**
 - JJ Gordon
 - J Gardener
 - J Ververs
 - PJ Keall
 - IB Mihaylov
 - W Li
 - I Kawrakow

VCU Virginia Commonwealth University