

Latest Advances in Digital Breast Tomosynthesis

Wei Zhao, Bo Zhao and Yue-houng Hu
Dept. of Radiology
State University of New York at Stony Brook

Limitation of Mammography


- False diagnoses due to superposition of breast tissue
- Solution: 3D imaging
 - Breast CT: Cone-beam CT (180°) with dedicated prone table
 - Breast tomosynthesis: modified mammo gantry with limited angular range (15°-60°)

Breast Tomosynthesis

- Fast clinical transition
- Limited angular range: 15~60°
- Slice thickness: 1mm
- Projections: 11- 50
- Dose: ~1 - 2 times single mammogram




Mammogram



Biopsy proven cancer

Tomosynthesis



Courtesy: B. Ren, Hologic


Different DBT Design

Company	System param.	Views	Detector	Scan time (s)	Recon.
GE	$\pm 20^\circ$, step/shoot (MGH, Sunnybrook)	15	CsI/a-Si 100 μm	15-23	MLEM SART
	$\pm 30^\circ$, step/shoot (U. of Michigan)	21		7 s	SART
Hologic	$\pm 7.5^\circ$, continuous Multiple sites	11	a-Se, 70 μm 2x2 binning	10	FBP
Siemens	$\pm 22^\circ$, step/shoot Duke, Malmö, SUNY	25	a-Se, 85 μm	12.5/20 Bin/full	FBP
Dexela	$\pm 12-20^\circ$ (U. of Virginia)	13	Fiber optic coupled CCD	30 s	MLEM
X-counter	$\pm 13^\circ$	48	Gas counting, 48 slit, 60 μm		FBP iterative
Sectra	$\pm 5.5^\circ$, 6-site trial	21	Si counting, 21 slit, 50 μm	3-8 s	


Factors Affecting Image Quality

- Detector performance
 - DQE at low dose: $1/N_{\text{View}}$
 - Temporal performance: lag and ghosting
- Imaging geometry
 - Angular range, number of views
 - Focal spot blur: continuous tube travel
- Reconstruction algorithm
 - Analytical: FBP (filter design)
 - Iterative





Detector Performance







Imaging Geometry: spatial vs. frequency domain




B. Zhao and W. Zhao, Med. Phys. 2008






MTF – Dependence on Filters




Impact of Angular Range on In-Plane MTF




• Increasing angular range improves MTF at low frequencies


In-Plane MTF Measurement



In-Plane MTF




Dependence on Angular Range



ACR phantom imaged at 28kVp, W/Rh, 1.7mGy
Recon: FBP with ST filter


Reconstruction filter designs

- **Ramp + limited angle:**
loss of low frequency (breast density)
- **Modified filters:**
recovers density, but increased out-of-plane blur
- **Iterative recon:**
improvement at the cost of computation



J. Zhou et. al. Med. Phys. 2007; Ludwig et. al. IWDM2008; B. Ren, et. al. SPIE 2009

Reconstruction filter comparison



FBP with Ramp





SIRT





Polynomial

Adapted from Ludwig et. al. (IWDM 2008)

Clinical Application of DBT

- Screening/diagnosis
- Contrast enhancement: angiogenesis
- Multimodality imaging:
 - DBT + automated ultrasound (AUS)
 - P. Carson et. al. U. Michigan
 - DBT + optical tomography
 - Boas et. al. MGH
 - DBT + limited angle SPECT
 - M. Williams et. al., U. Virginia






Further Development for DBT Systems

- **Gantry:**
 - Faster scan, more views
- **Generator/tube:**
 - dual-energy, shorter pulse/less focal spot blur
- **Detector:**
 - Better low dose performance, faster readout
 - photon-counting

New detectors: Improve low dose performance



Summary

- Physics of DBT: 2D→3D
 - Optimization of acquisition and recon.
- New clinical applications:
 - Contrast enhanced CEDBT
 - Multimodality imaging
- Challenges for further system development

Acknowledgements

- Financial support:
 - Siemens AG
 - Army Breast Cancer Research Program
 - NIH
- Contributing slides:
 - Drs. Baorui Ren, James Mainprize, Mitch Goodsitt, Ann-Katherine Carton, Mats Danielsson