



Advanced Ultrasound Imaging in Interventional Medicine

Emad M. Boctor, Ph.D.

Assistant Professor of Radiology, Division of Medical Imaging Physics,
The Russell H. Morgan Department of Radiology
and Radiological Science.

eboctor@jhmi.edu

AAPM 2010




Disclaimer

- Overview and understanding, not comparison
- Not every possible work will be discussed (lack of time)
- Not every group or individual will be covered (lack of time)
- There is no financial interest with the companies mentioned in this presentation

Driving Application: Liver Ablation



- Hepatocellular Carcinoma (HCC): 1M cases per year (worldwide)
- The most frequent hepatic malignancy in USA is metastatic disease from colorectal cancer
- Resection -- 5 year survival rates between 25% and 55%
- Most patients do not qualify for resection
- un-resectable liver tumors are ablated under ultrasound guidance




Thermal Ablation of Liver Tumors

Problems with the Free-Hand Approach


- Dependent on physician accuracy
- Often requires multiple passes
- Unsuccessful ablation rate = 5%
- Inconsistent
- Not repeatable
- Post-ablation evaluation




Situations in Which 3DUS Guidance May Be Most Useful




Proposed Solutions



A Dual-Armed Robotic System



• Boctor E.M., Fischer G., Choti M.A., Fichtinger G., and Taylor R.H. "Dual-Armed Robotic System for Intraoperative Ultrasound Guided Hepatic Ablative Therapy: A Prospective Study", IEEE 2004 International Conference on Robotics and Automation, in proceedings, pp. 377-382.



Ultrasound Calibration

$$\begin{pmatrix} 0 \\ 0 \\ 0 \\ 1 \end{pmatrix} = c_{T_x}^T T_B^B R_{T_x} \begin{pmatrix} x_{w1} \\ x_{w2} \\ x_{w3} \\ 1 \end{pmatrix}$$

Courtesy of R. Prager

- Prager et al. "Rapid calibration for 3-D freehand ultrasound" UMB 1998.
- Galloway et al, UMB 2001, Abolmaesumi et al, MICCAI 2006, Khamene et al, MICCAI 2005, ...

Closed form formulation

$$B = B_2^{-1} B_1$$

$$A = A_2 A_1^{-1}$$

↓

$$AX = XB$$

- Bector E.M., Viswanathan A., Choti M.A., Taylor R.H., Fichtinger G., and Hager G.D. "A Novel Closed Form Solution For Ultrasound Calibration", ISBI 2004, in proceedings, pp 527-530.
- Danilidis et al., IJRR 1999.

How to solve $AX=XB$?

$$R_a R_b = R_c R_d$$

$$R_a t_x + \lambda t_a = R_c t_b + t_x$$

$$(R_a \otimes R_b) \text{vec}(R_x) = \text{vec}(R_x)$$

$$(I_3 \otimes t'_b) \text{vec}(R_x) + (I_3 - R_a) t_x - \lambda t_a = 0$$

$$\begin{bmatrix} I_3 - R_a \otimes R_b & 0_{9 \times 3} & 0_{9 \times 3} \\ I_3 \otimes t'_b & I_3 - R_a & -t_a \end{bmatrix} \begin{pmatrix} \text{vec}(R_x) \\ t_x \\ \lambda \end{pmatrix} = \begin{pmatrix} 0_{9 \times 1} \\ 0_{3 \times 1} \end{pmatrix}$$

$$UV + VW = T$$

$$(U \otimes I + I \otimes W) \text{vec}(V) = \text{vec}(T)$$

$$\text{vec}(CDE) = (C \otimes E^T) \text{vec}(D)$$


Patient Specific in vivo Calibration

If only we could estimate "A" without phantom...


$$AX = XB$$

- Bector et al., MICCAI 2005, and SPIE 2006.
- Wein and Khamene SPIE 2008.
- Barratt et al., MICCAI 2005.

Ablation under US Guidance is Blind




B-mode image



Gross-pathology


Elastography (Pioneered by Ophir, Sarvazyan, Bamber, Varghese, Hall, Emelianov, ...)



2D representation of strain based imaging model. Before compression: the overlay represents 1D cascaded particles with uniform spacing. After compression: the overlay represents two groups of particle spacing. Small spacing (light green) indicating soft tissues moved more (high strain) than the hard tissue (low strain).


Stress-Strain Measurements

- Elasticity changes are immediate and permanent
- Cooked and raw liver can always be told apart
 - Young's modulus ratio is ~10
 - Stress is linear below ~5% strain





Ex vivo Imaging Study

- Supporting gelatin
- Fiducial markers in transparent gelatin
- Radionics single-rod ablator device
- Ellipsoidal ablation along the needle shaft
- Large ablation in short time by using cool-tip technology




Registration between Elasticity Image and Gross-pathology


- Supporting gelatin
- Fiducial markers in transparent gelatin
- Radionics single-rod ablator device
- Ellipsoidal ablation along the needle shaft
- Large ablation in short time by using cool-tip technology



The Liver Samples







Strain Results



B-mode image shows ex-vivo liver boundaries embedded in gel based medium. It is not possible to differentiate the ablated area from B-mode. Strain is generated from differentiating a displacement map in the axial direction. Strain provides clear evidence of the presence of hard lesion, which is in agreement with the gross pathology picture.

Serial Segmentation Pipeline





- ### Moving from ex vivo to in vivo
- Real-time strain imaging or rapid interactive rate
 - Robustness to uncontrolled motion
 - High resolution, SNR and CNR
 - 2D (or 3D) extension
 - High axial compression
 - Insensitivity to signal decorrelation

Dynamic programming approach

Amplitude similarity

$$\Delta(i, d) = |g(i) - g'(i + d)|$$


Smoothness

$$S(d_i, d_{i-1}) = (d_i - d_{i-1})^k$$

Recursive cost function

$$C(i, d_i) = \min_{d_{i-1}} \{C(i-1, d_{i-1}) + wS(d_i, d_{i-1})\} + \Delta(i, d_i)$$

$$M(i, d_i) = \arg \min_{d_{i-1}} \{C(i-1, d_{i-1}) + wS(d_i, d_{i-1})\}$$






• Hager et al., PAMI 2003; Hall et al., US IEEE Sym. 2006;

Contrast to Noise Ratio


$$CNR = \frac{C}{N} = \sqrt{\frac{2(\bar{s}_b - \bar{s}_i)^2}{\sigma_b^2 + \sigma_i^2}}$$

- Target window is fixed on the lesion
- Background window, is moved across the strain image




• Rivaz et al., TMI 2008.


Dynamic Programming Elastography vs. Normalized Cross-correlation Methods




2D Dynamic Programming Elastography



Freehand Palpation of Resected Prostate



In vivo Patient Studies




Rivaz et al., MICCAI 2008

Challenges and Possible Solutions


- From 2D to 3D displacement
- Effective and rapid visualization
- **Optimal** real-time elasticity imaging
- Large animal model for in vivo validation

3D Elasticity Imaging of Ablation



Rivaz et al., MICCAI 2008

Elasticity-based Volume Rendering of 3DUS B-mode data



Ray Casting Volume Rendering Pipeline Based on Strain Data as Opacity Volume

Elasticity-based Volume Rendering of 3DUS B-mode data



Acknowledgements

Research Collaborators:

Gregory Hager (Research Director)
 Russell Taylor (Director)
 Iulian Iordachita (CS, Hopkins)
 Gabor Fichtinger (CS, Queens)
 Roger Ghanem (CE, USC)
 Gregory Chirikjian (ME, Hopkins)
 Jin Kang (ECE, Hopkins)
 Jim Spicer (ME, Hopkins)
 Colleagues and Staff in the ERC and LCSR

Funding:

- SIEMENS Corporate Research, pre-doctoral fellowship
- NSF EEC 9731478, ERC Center Grant
- BCRF

Clinical Collaborators:

- Michael Choti, M. Awad, L. Assumcao, M.DeOliviera (Surgery, Hopkins)
- Ted DeWeese, Danny Song, R. Zellars (Rad-On, Hopkins)
- Ron Kikinis (Surgical Planning Lab, BWH)

Industrial Collaborators:

- Frank Sauer, Ali Khamene, Wolfgang Wein (SIEMENS Corporate Research)
- Shelby Brunke, SIEMENS Ultrasound Dept.
- CMS, Burdette Medical Systems
- Intuitive Surgical Inc.
- Aloka Ultrasound



Thank you !



Follow-up:
Emad Bactor
ebactor@jhmi.edu

