

Current Status of HIFU Therapy for Treatment of Benign and Malignant Tumors of the Abdomen, Pelvis and Bone

Joo Ha Hwang, MD, PhD

Department of Medicine, Division of Gastroenterology
 Department of Radiology
 Applied Physics Laboratory, Center for Industrial and Medical Ultrasound
 University of Washington, Seattle, USA

UW Medicine
 SCHOOL OF MEDICINE



Objectives

- Overview of current HIFU clinical systems
- Overview of current clinical applications
 - Uterine fibroids
 - Pancreatic tumors
 - Liver tumors
 - Renal cell carcinoma
 - Bone metastases
 - Breast cancer
 - Thyroid/parathyroid tumors
- Discuss future clinical applications

Clinical HIFU Systems US-guided



YDME, Beijing



HIAFU, Chongqing




Shanghi A&S



Thereclion, France

Ultrasound targeting

- Ultrasound-guided HIFU devices do not provide monitoring of lesion development (other than hyperecho from boiling)
- Methods for estimating *in situ* intensity exist
- Methods for monitoring HIFU therapy are in development
 - Thermometry
 - ARFI
 - Elastography



Clinical HIFU Systems MR-guided



InSightec, Israel




Philips, US/France



Temperature mapping

Proton Resonance Frequency shift

Phase Imaging



$$\Delta T = \frac{\Delta\Phi}{\alpha \cdot \gamma \cdot T_E \cdot B_0}$$

$\gamma = 2\pi \cdot 42.56 \text{ MHz/T}$
 $\alpha = 0.0101 \text{ ppm/}^\circ\text{C}$
 $T_E \sim 20\text{ms}$
 $B_0 = 1.5 \text{ T}$

Gyromagnetic Ratio
 Water Frequency Shift
 Echo Time
 Magnetic Field

Courtesy of Philips

Real Time Feedback

Reliable necrosis volume



Thermal map & Dose map

Real time visualization
+ Feedback

$T > 57^\circ\text{C}^*$ or Dose $> 240 \text{ EM}$

Stop heating

Reliable necrosis volume
No a-priori knowledge needed
Simple and robust



Non-perfused volume

* Applies to the border of the cell. Temperatures at the center are higher, especially for larger cells.

Courtesy of Philips

Clinical Applications for HIFU Ablation

- Indications: Non-invasive ablation of solid benign or malignant tumors
- Requirement:
 - Acoustic window – critical
- Non-oncologic applications
 - Uterine fibroids

Treatment of Uterine Fibroids with HIFU

Uterine Fibroids



Pre-HIFU

Post-HIFU

Courtesy of Dr. Hu

HIFU of Fibroids - Fertility

- Uterine leiomyomas (fibroids) affect ~25% of women of reproductive age
- Large nonhysteroscopically resectable submucosal and intramural fibroids can cause cavity distortion impacting fertility
- MRgFUS
 - 54 pregnancies in 51 women
 - Live births (41%)
 - Spontaneous abortion (28%)
 - Elective abortion (11%)
 - Ongoing pregnancies beyond 20 wks (20%)
 - Vaginal delivery rate – 64%

Hanstede et al. Fertil Steril 2007
Rabinovici et al. Fertil Steril 2008

Oncologic Applications for HIFU Ablation

- Oncologic indications:
 - Palliation
 - Local tumor control
 - Poor surgical candidate
 - Patient refuses surgery
- Oncologic applications:
 - Pancreatic cancer
 - Liver tumors (should be below costal margin)
 - Renal cell carcinoma
 - Osteosarcoma/soft tissue sarcomas/bone metastasis
 - Breast cancer

HIFU Ablation for Palliation of Advanced Pancreatic Cancer


Pancreatic Cancer

- 4th leading cause of cancer deaths in the US
- >42,000 new cases – >35,000 deaths in 2010
- 80-90% are “unresectable” when diagnosed
- Poor outcomes
 - Median survival without therapy: 3 months
 - Median survival with therapy: 6-12 months
- Palliation of symptoms is important
 - Pain relief

Treatment of Pancreatic Cancer



Pancreatic Cancer



HIFU Ablation of Liver Tumors

HIFU Treatment of Liver Tumors

- Acoustic window is a problem
- Treatment through the ribs has been reported
- Treatment after rib resection
- Randomized study (Wu et al.)
 - TACE vs. HIFU+TACE (50 patients)
 - HIFU treatment 2-4 weeks following TACE

outcome	TACE	TACE+HIFU	p-value
Median survival	4.0 months	11.3 months	0.004
6-month survival	13.2%	80.5%	0.002
1-year survival	0%	42.9%	<0.001

Liver Tumors



Pre-HIFU




3 months post- HIFU

Courtesy of Dr. Hu

HIFU Ablation of Renal Cell Carcinoma

Renal Cell Carcinoma

- Small renal tumors are being more frequently identified
- Surgical resection has significant morbidity
- Many pts are poor surgical candidates
- Other less invasive procedures are needed



Marberger et al. BJU 2005

HIFU Ablation of Bone Metastasis

Bone Metastasis

- Bone is a common site for metastatic disease
 - Prostate CA and Breast CA
- Pain from bone metastasis is the most common cause of cancer pain
- Current treatment options:
 - Analgesics
 - Chemotherapy
 - Biphosphonates
 - Local therapies:
 - Radiation (no relief in 20-30%)
 - Surgery
 - RFA




Feasibility Study

Patient population:

36 treatments in 31 patients were conducted, targeting 32 metastatic lesions


Patient tumor characteristics:

- Treated bone mets from **primary tumor type**: renal, colorectal, lung, breast, prostate and other cancers
- Treated **lesion locations**: iliac bone, ischium bone, sacrum, femur, scapula, humerus, clavicle
- Treated **lesion type**: both osteolytic and osteoblastic



Lieberman B et al. Pain palliation in patients with bone metastases using MR-guided focused ultrasound surgery: a multicenter study. Ann Surg Oncol. 2009 Jan;16(1):140-6. Epub 2008 Nov 11.

Bone Metastasis – Palliation of Pain




Catane et al. Ann Oncol 2007

HIFU Ablation of Breast Cancer

Breast Cancer – Phase II Study Results.


- Evaluation of Safety and Efficacy
- Total 195 patients treated, 30 were treated with contrast enhanced planning images



The contrast enhanced subcutaneous images

Mean necrosis of ablated tumor = $96.9 \pm 4\%$

- 15 patients had 100% necrosis
- 3 patients had < 95% necrosis





Post-treatment pathologic results

Furusawa H, Namba K, Thomsen S, Miyama F, Berndt A, Tanaka G, Yasuda Y, Hatahara H. Magnetic Resonance-Guided Focused Ultrasound Surgery of Breast Cancer: Reliability and Effectiveness. J Am Coll Surg. 2006; 203(1):54-63

Courtesy of Breastopia Namba Medical Center, Miyazaki, Japan


HIFU Ablation of Thyroid/Parathyroid Tumors

Clinical studies



Courtesy of Theracision

Gross Specimen



Courtesy of Theracision

HIFU: Future

The Future of HIFU

- Expanded indications for therapy
- Enhanced tumor-specific immunity
- Drug delivery
 - Targeted delivery (Heat activation)
 - Enhancement of vascular permeability
 - Enhanced penetration
- Improved treatment monitoring
 - Thermometry
 - Elastography
 - Radiation Force Imaging
- Device approvals

Summary

- Overview of current clinical systems
- Overview of current clinical applications
 - Uterine fibroids
 - Pancreatic tumors
 - Liver tumors
 - Renal cell carcinoma
 - Bone metastases
 - Breast cancer
 - Thyroid/parathyroid tumors
- Discussed future clinical applications

Thank You